

User Manual

PAX-327-C8 Series

27" 4K True-Flat Medical Monitor

ADVANTECH

Enabling an Intelligent Planet

Copyright

The documentation and the software included with this product are copyrighted 2021 by Advantech Co., Ltd. All rights are reserved. Advantech Co., Ltd. reserves the right to make improvements in the products described in this manual at any time without notice. No part of this manual may be reproduced, copied, translated, or transmitted in any form or by any means without the prior written permission of Advantech Co., Ltd. The information provided in this manual is intended to be accurate and reliable. However, Advantech Co., Ltd. assumes no responsibility for its use, nor for any infringements of the rights of third parties that may result from its use.

Trademarks

All trademarks and registered trademarks are property of their respective owners.

Disclaimer

Although every attempt has been made to achieve technical accuracy in this manual, we assume no responsibility for errors that may be found. Our goal is to provide you with the most accurate and usable documentation possible. To assist us with improving this manual, we welcome all comments and constructive criticism. Please send all feedback in writing to support@advantech.com.

Important

Read this user manual carefully to familiarize yourself with the correct, safe, and effective usage procedures. Additionally, we recommend that you retain this manual for future reference.

Importante

Lea detenidamente este manual del usuario para familiarizarse con los procedimientos de uso eficaces y seguros. Conserve este manual para futuras consultas.

重要

ご使用前には必ず取扱説明書をよくお読みになり、正しくお使いください。この取扱説明書は大切に保管してください。

User Instructions

This document combines text and illustrations to provide a comprehensive overview of the system. The information is presented as a sequential step of actions, allowing the user to learn directly how to use the device.

The text provides explanations and instructs the user step-by-step in the practical use of the product with short and clear instructions in an easy-to-follow sequence.

Additional Information and Assistance

Contact your distributor, sales representative, or an Advantech customer service center for technical support if you need additional assistance. Please have the following information ready before calling:

- Product name and serial number
- Description of your peripheral attachments
- Description of your software (operating system, version, application software, etc.)
- A complete description of the problem
- The exact wording of any error messages

Note: This equipment is a source of electromagnetic waves. Before use, ensure that no EMI-sensitive devices that may malfunction are in the surrounding area.

Manufacturer

Advantech Co., Ltd.

No.1, Alley 20, Lane 26, Rueiguang Road Neihu District, Taipei,
Taiwan 114, R.O.C.

Tel: (02) 2792-7818

Visit the Advantech website at www.advantech.com or www.advantech.com.tw for additional information.

Contents

Chapter 1	Introduction.....	1
1.1	Introduction	2
1.2	Features	2
1.3	Specifications	3
1.3.1	Display	3
1.3.2	Input / Output Terminals	5
1.3.3	General	8
1.3.4	Environmental condition.....	9
1.3.5	Optional Advantech Video Archiving Solution (AVAS).....	9
1.3.6	I/O Connectors.....	10
Chapter 2	Installation.....	13
2.1	Unpacking	14
2.2	System Controls	14
	Figure 2.1 PAX-327-C8 Front View	14
2.2.1	I/O Ports.....	15
2.2.2	Display	17
2.2.3	OSD Key Functions	18
2.2.4	Input/Output Signal Ports.....	18
2.2.5	Interface on AVAS Back Panel	21
2.3	Viewing Angle	21
2.4	Positioning.....	22
2.5	Connecting AC/DC Power.....	22
2.6	Connecting Video.....	22
Chapter 3	OSD Menu	27
3.1	OSD Menu Controls	28
3.1.1	Display Mode	28
3.1.2	Display Function	29
3.1.3	Select Region.....	30
3.1.4	Picture.....	31
3.1.5	Color	32
3.1.6	Advance	33
3.1.7	Input.....	34
3.1.8	Audio (Option).....	35
3.1.9	Other	36
3.1.10	Information	37
3.1.11	OSD Input Source.....	37
3.1.12	PIP Table	37
Chapter 4	Safety Information	39
4.1	Symbols	40
4.2	FCC Compliance	40
4.2.1	Class A.....	40
4.3	CE Marking	41
4.4	VCCI (Voluntary Control Council for Interference)	41
4.5	Disposal of Waste.....	41
4.6	Safety and Maintenance	42
4.6.1	Radio Frequency.....	43

4.6.2	Keep Away From Windows.....	43
4.6.3	Storage Temperature.....	43
4.6.4	Image Persistence.....	43
4.6.5	Placement and Adjustment of the Monitor.....	44

Appendix A **Cleaning..... 45**

A.1	Cleaning the LCD Panel	46
A.2	Cleaning the Cabinet	46

Appendix B **Troubleshooting..... 47**

B.1	Non-Responsive Power Button.....	48
B.2	No Picture	48
B.3	No Video	48
B.4	Image Persistence	48
B.5	Unstable, Unfocused, or Swimming Image.....	48
B.6	Incorrect Display Image Size	48

Chapter 1

Introduction

1.1 Introduction

PAX-327-C8 is a 27-inch TFT LCD monitor designed for medical imaging display applications. With the digital imaging and communication in medicine (DICOM) gray-scale standard display function (GSDF) and 14-bit gamma correction for white balance, PAX-327-C8 is optimized for displaying high-quality medical images and videos.

1.2 Features

- True-flat design with AR filter (optional)
- High brightness and high contrast ratio
- Integral LED backlight
- Wide viewing angle
- 27-inch wide-screen with UHD resolution
- HDMI 2.0 supports up to 3840 x 2160 @ 60 Hz resolution
- DisplayPort 1.2 supports up to 3840 x 2160 @ 60 Hz resolution
- SDI supports up to 3840 x 2160 @ 60 Hz resolution (optional)
- DVI supports up to 3840 x 2160 @ 30 Hz resolution
- Dual modality connectivity
- Automatic detection of input source for signal switching
- Video enhancement features include
 - 14-bit LUT processing
 - DICOM Part 14 GSDF compliance
 - DICOM clear and blue modes
- Variable gamma and DICOM (color temperature) modes
- Variable picture control
- Variable image size
- Variable multi-picture display mode
- Built-in Advantech AVAS video solutions (optional)
- Integrated touchscreen (optional)
- Wide input connectivity for maximum flexibility

1.3 Specifications

1.3.1 Display

Model Name	Specifications			
	PAX-327-C8A-NAS PAX-327-C8A-NAR	PAX-327-C8E-NAS PAX-327-C8E-NAR	PAX-327-C8D-NAS PAX-327-C8D-NAR	PAX-327-C8F-NAS
Screen Size (Active Area)	27" (596.74(H) x 335.66(V) [mm])			
Aspect Ratio	16:9			
Number of Pixels	3840 (H) x 2160 (V)			
Pixel Pitch	0.1554(H) x 0.1554(V) [mm]			
Displayable Colors	1.07 billion (10-bit)			
Brightness (Typical)	350 cd/m ²			
Contrast Ratio (Typical)	1300:1			
Display Mode	IPS			
Response Time (Typical)	14ms (Gray to Gray)			
Viewing Angle (H/V)	178/178			
Touchscreen (Optional)	Projected capacitive			

Model Name	Specifications			
	PAX-327-C8A-HAR-AI	PAX-327-C8A-HAS PAX-327-C8A-HAR	PAX-327-C8E-HAS PAX-327-C8E-HAR	PAX-327-C8D-HAS PAX-327-C8D-HAR
Screen Size (Active Area)	27" (596.16(H) x 335.34(V) [mm])			
Aspect Ratio	16:9			
Number of Pixels	3840 (H) x 2160 (V)			
Pixel Pitch	0.15525(H) x 0.15525(V) [mm]			
Displayable Colors	1.07 billion (10 bit)			
Brightness (Typical)	800 cd/m ²			
Contrast Ratio (Typical)	1,000:1			
Display Mode	AHVA			
Response Time (Typical)	16 ms (Tr+Tf)			
Viewing Angle (H/V)	178/178			

Touchscreen (Optional)	Projected capacitive
-------------------------------	----------------------

Model Name	Specifications		
	PAX-327-C8L-HAS PAX-327-C8L-HAR	PAX-327-C8F-HAS	PAX-327-C8J-HAS
Screen Size (Active Area)	27" (596.16(H) x 335.34(V) [mm])		
Aspect Ratio	16:9		
Number of Pixels	3840 (H) x 2160 (V)		
Pixel Pitch	0.15525(H) x 0.15525(V) [mm]		
Displayable Colors	1.07 billion (10 bit)		
Brightness (Typical)	800 cd/m ²		
Contrast Ratio (Typical)	1,000:1		
Display Mode	AHVA		
Response Time (Typical)	16 ms (Tr+Tf)		
Viewing Angle (H/V)	178/178		
Touchscreen (Optional)	Projected capacitive		

1.3.2 Input / Output Terminals

Item		Specifications			
		PAX-327-C8A-NAS PAX-327-C8A-NAR	PAX-327-C8E-NAS PAX-327-C8E-NAR	PAX-327-C8D-NAS PAX-327-C8D-NAR	PAX-327-C8F-NAS
DC Power Input	DC Jack x 1	DC24V, 4.16A(MAX) or 5A(MAX)			
DC Power Output	DC Jack x 1	DC5V or DC12V			
RS-232C	D-Sub Jack x 1	Service Port			
DVI Input	DVI Jack x 1 (Type D)	Digital RGB: TMDS MAX: 3840 x 2160@30Hz			
DVI Output	DVI Jack x 1 (Type D)	Digital RGB: TMDS MAX: 3840 x 2160@30Hz			
HDMI 2.0 Input	HDMI Jack x 1 (Type A)	Digital RGB: TMDS MAX: 3840 x 2160@60Hz			
HDMI 2.0 (For Ndcoder)	HDMI Jack x 1 (Type A)	Digital RGB: TMDS MAX: 3840 x 2160@60Hz			
HDMI 2.0 Output	HDMI Jack x 1 (Type A)	Digital RGB: TMDS MAX: 3840 x 2160@60Hz			
DP 1.2 Input	DP Jack x 1	DisplayPort 1.2 MAX: 3840 x 2160@60Hz			
DP 1.2 Out-put	DP Jack x 1	-			
USB Up	USB-B Jack x 1	USB 2.0 (For Touch Screen)			
SDI Quad 3G Input	BNC Jack x 4	-	-	Serial Input Common Mode: 2.2V[p-p](75Ω) MAX: 3840 x 2160@60Hz	
SDI Quad 3G Output	BNC Jack x 4	-	-	Serial Output Common Mode: 2.2V[p-p](75Ω) MAX: 3840 x 2160@60Hz	
SDI 12G Input	BNC Jack x 1	-	Serial Input Com- mon Mode: 2.2V[p-p](75Ω) MAX: 3840 x 2160@60Hz	-	-
SDI 12G Output	BNC Jack x 1	-	Serial Output Com- mon Mode: 2.2V[p-p](75Ω) MAX: 3840 x 2160@60Hz	-	-
DVI Output (For Display)	DVI Jack x 1 (Type D) (Internal Connection)	-	-	-	Digital RGB: TMDS MAX: 1920 x 1080@60Hz
LAN	RJ45 Jack x 1	-	-	-	1GbE
SFP+	Module	-	-	-	SFP+ 10Gbps

Item		Specifications			
		PAX-327-C8A-HAR-AI	PAX-327-C8A-HAS PAX-327-C8A-HAR	PAX-327-C8E-HAS PAX-327-C8E-HAR	PAX-327-C8D-HAS PAX-327-C8D-HAR
DC Power Input	DC Jack x 1	DC24V, 3.75A(MAX)	DC24V, 6.25A(MAX) or 5A(MAX)		
DC Power Output	DC Jack x 1	DC5V or DC12V			
RS-232C	D-Sub Jack x 1	Service Port			
DVI Input	DVI Jack x 1 (Type D)	Digital RGB: TMDS MAX: 3840 x 2160@30Hz			
DVI Output	DVI Jack x 1 (Type D)	Digital RGB: TMDS MAX: 3840 x 2160@30Hz			
HDMI 2.0 Input	HDMI Jack x 1 (Type A)	Digital RGB: TMDS MAX: 3840 x 2160@60Hz			
HDMI 2.0 (For Ndcoder)	HDMI Jack x 1 (Type A)	Digital RGB: TMDS MAX: 3840 x 2160@60Hz			
HDMI 2.0 Output	HDMI Jack x 1 (Type A)	Digital RGB: TMDS MAX: 3840 x 2160@60Hz			
DP 1.2 Input	DP Jack x 1	DisplayPort 1.2 MAX: 3840 x 2160@60Hz			
DP 1.2 Out-put	DP Jack x 1	DP(1.2) MST or Clone			
USB Up	USB-B Jack x 1	USB 2.0 (For Touch Screen)			
SDI Quad 3G Input	BNC Jack x 4	-	-	-	Serial Input Com- mon Mode: 2.2V[p-p](75Ω) MAX: 3840 x 2160@60Hz
SDI Quad 3G Output	BNC Jack x 4	-	-	-	Serial Output Com- mon Mode: 2.2V[p-p](75Ω) MAX: 3840 x 2160@60Hz
SDI 12G Input	BNC Jack x 1	-	-	Serial Input Com- mon Mode: 2.2V[p-p](75Ω) MAX: 3840 x 2160@60Hz	-
SDI 12G Output	BNC Jack x 1	-	-	Serial Output Com- mon Mode: 2.2V[p-p](75Ω) MAX: 3840 x 2160@60Hz	-
DVI Output (For Display)	DVI Jack x 1 (Type D) (Internal Connection)	-	-	-	-
LAN	RJ45 Jack x 1	-	-	-	-
SFP+	Module	-	-	-	-

Item		Specifications		
		PAX-327-C8L-HAS PAX-327-C8L-HAR	PAX-327-C8F-HAS	PAX-327-C8J-HAS
DC Power Input	DC Jack x 1	DC24V, 6.25A(MAX) or 5A(MAX)		
DC Power Output	DC Jack x 1	DC5V or DC12V		
RS-232C	D-Sub Jack x 1	Service Port		
DVI Input	DVI Jack x 1 (Type D)	Digital RGB: TMDS MAX: 3840 x 2160@30Hz		
DVI Output	DVI Jack x 1 (Type D)	Digital RGB: TMDS MAX: 3840 x 2160@30Hz		
HDMI 2.0 Input	HDMI Jack x 1 (Type A)	Digital RGB: TMDS MAX: 3840 x 2160@60Hz		
HDMI 2.0 (For Ncdecoder)	HDMI Jack x 1 (Type A)	Digital RGB: TMDS MAX: 3840 x 2160@60Hz		
HDMI 2.0 Output	HDMI Jack x 1 (Type A)	Digital RGB: TMDS MAX: 3840 x 2160@60Hz		
DP 1.2 Input	DP Jack x 1	DisplayPort 1.2 MAX: 3840 x 2160@60Hz		
DP 1.2 Output	DP Jack x 1	DP(1.2) MST or Clone		
USB Up	USB-B Jack x 1	USB 2.0 (For Touch Screen)		
SDI Quad 3G Input	BNC Jack x 4	Serial Input Common Mode: 2.2V[p-p](75Ω) MAX: 3840 x 2160@60Hz	Serial Input Common Mode: 2.2V[p-p](75Ω) MAX: 3840 x 2160@60Hz	-
SDI Quad 3G Output	BNC Jack x 4	Serial Output Common Mode: 2.2V[p-p](75Ω) MAX: 3840 x 2160@60Hz	Serial Output Common Mode: 2.2V[p-p](75Ω) MAX: 3840 x 2160@60Hz	-
SDI 12G Input	BNC Jack x 1	Serial Input Common Mode: 2.2V[p-p](75Ω) MAX: 3840 x 2160@60Hz	-	-
SDI 12G Output	BNC Jack x 1	Serial Output Common Mode: 2.2V[p-p](75Ω) MAX: 3840 x 2160@60Hz	-	-
DVI Output (For Display)	DVI Jack x 1 (Type D) (Internal Connection)	-	Digital RGB: TMDS MAX: 1920 x 1080@60Hz	
LAN	RJ45 Jack x 1	-	1GbE	
SFP+	Module	-	SFP+ 10Gbps	

Caution! Do not touch any of the signal inputs, signal outputs, or other connectors on the device while touching a patient.



External equipment intended for connection with the signal input, signal output, or other connectors must comply with the relevant IEC standard (e.g., IEC 60950 for IT equipment and IEC 60601-1 for medical electrical equipment). Additionally, all combination systems must comply with the IEC 60601-1 and IEC 60601-1-1 specifications.

If in doubt, contact a qualified technician or your local sales rep.

Note! Although Advantech products are subject to change without notice, the following specifications were accurate at the time of manual publication:



1. Type of protection against electric shock: Class I equipment
2. Degree of protection against electric shock: Not classified - no applied parts
3. Classification according to the degree of protection against water ingress: IPX0, ordinary equipment
4. This equipment is not suitable for use in the presence of flammable anesthetics or oxygen
5. Mode of operation: continuous operation

1.3.3 General

Model Name	Specifications			
	PAX-327-C8A-NAS PAX-327-C8A-NAR	PAX-327-C8E-NAS PAX-327-C8E-NAR	PAX-327-C8D-NAS PAX-327-C8D-NAR	PAX-327-C8F-NAS
Power Supply	DC24V, 4.16A(MAX) or 5A(MAX)			
Power Consumption (Typical)	39.8 W	64.3 W	67.9 W	72.2 W
Dimensions (Without Stand) [W mm x H mm x D mm]	663.83 x 418.11 x 73.55 [mm]			
Weight (Without Stand)	9.32 kg	9.52 kg	9.78 kg	10.16 kg

Model Name	Specifications			
	PAX-327-C8A-HAR-AI	PAX-327-C8A-HAS PAX-327-C8A-HAR	PAX-327-C8E-HAS PAX-327-C8E-HAR	PAX-327-C8D-HAS PAX-327-C8D-HAR
Power Supply	DC24V, 3.75 A max.	DC24V, 6.25 A or 5 A max.		
Power Consumption (Typical)	66.8 W	66.8 W	91.8 W	74.1 W
Dimensions (Without Stand) [W mm x H mm x D mm]	663.83 x 418.11 x 73.55 [mm]			
Weight (Without Stand)	9.52 kg	9.52 kg	9.72 kg	9.98 kg

Model Name	Specifications		
	PAX-327-C8L-HAS PAX-327-C8L-HAR	PAX-327-C8F-HAS	PAX-327-C8J-HAS
Power Supply	DC24V, 6.25 A or 5 A max.		
Power Consumption (Typical)	84.9 W	100.4 W	92.6 W
Dimensions (Without Stand) [W mm x H mm x D mm]	663.83 x 418.11 x 73.55 [mm]		
Weight (Without Stand)	9.99 kg	10.36 kg	9.90 kg

1.3.4 Environmental condition

		Model Name
		PAX-327-C8 series
Temperature	Operating	0 ~ 40 °C
	Storage	-20 ~ 60 °C
Humidity	Operating	Ta = 40 °C/104 °F, 90% RH (non-condensing)
	Storage	5 ~ 90%
Pressure	Operating	500 ~ 1013 hPa
	Storage	500 ~ 1013 hPa

Note! * The specifications are subject to change without notice



1.3.5 Optional Advantech Video Archiving Solution (AVAS)

LAN	1 x RJ45 jack	1 x GbE
-----	---------------	---------

SFP+	Module	SFP+ 10 Gbps
------	--------	--------------

Built-in LEDs indicate the AVAS status, as explained below.


1. Blinking green: operating as normal
2. Blinking orange: warning notification
3. Continuous orange: FPGA is starting up
4. Blinking red: error notification
5. Off: no power supplied to the AVAS

Note! *Disconnecting the power supply cable will terminate power to the AVAS modules. The equipment should be installed near a power outlet that is easily accessible.*




1.3.6 I/O Connectors



PAX-327-C8A-NAS, PAX-327-C8A-NAR




PAX-327-C8E-NAS, PAX-327-C8E-NAR




PAX-327-C8D-NAS, PAX-327-C8D-NAR




PAX-327-C8F-NAS




PAX-327-C8A-HAS, PAX-327-C8A-HAR, PAX-327-C8A-HAR-AI




PAX-327-C8E-HAS, PAX-327-C8E-HAR




PAX-327-C8D-HAS, PAX-327-C8D-HAR



PAX-327-C8L-HAS, PAX-327-C8L-HAR




PAX-327-C8F-HAS



Note! For connecting the SDI cable, the 3G-SDI 1-lane input for HD/3G signals should be connected to Input Port 1. Meanwhile, the 3G-SDI 4-lane input for UHD signals should be connected to Input Ports 1, 2, 3 and 4.



PAX-327-C8J-HAS



Chapter 2

Installation

2.1 Unpacking

Before unpacking the monitor, prepare a suitable, level, and clean workspace near a power outlet. The monitor should be installed in a location with sufficient airflow and away from direct sunlight. After unpacking, check to ensure that the following items were included in the shipment:

- 1 x PAX-327 LCD monitor
- 1 x PAX-327 user manual

If any of the above items are missing or damaged, contact your dealer or sales representative immediately. Retain the shipping carton and packing material for storing or transporting the monitor in the future.

2.2 System Controls

The PAX-327 series monitors are designed to offer easy and convenient access to all control keys and peripheral ports. Before installation, take the time to familiarize yourself with the system controls and I/O ports.





Figure 2.1 PAX-327-C8 Front View




2.2.1 I/O Ports


PAX-327-C8A-NAS, PAX-327-C8A-NAR




PAX-327-C8E-NAS, PAX-327-C8E-NAR




PAX-327-C8D-NAS, PAX-327-C8D-NAR




PAX-327-C8F-NAS





PAX-327-C8A-HAS, PAX-327-C8A-HAR, PAX-327-C8A-HAR-AI



PAX-327-C8E-HAS, PAX-327-C8E-HAR




PAX-327-C8D-HAS, PAX-327-C8D-HAR




(9)


PAX-327-C8L-HAS, PAX-327-C8L-HAR



PAX-327-C8F-HAS



PAX-327-C8J-HAS



2.2.2 Display

All PAX-327-C8 series monitors feature a 27-inch TFT LCD panel that supports a max. resolution of 3840 x 2160 @ 60Hz.

2.2.3 OSD Key Functions

Key Name	Function
Power Key	
	Power On/Off
Menu/Enter Key	Select main OSD menu Select sub-menu Apply adjustments (save) Wake up in DPMS mode
	
Exit/Input Key	Exit/Return to previous menu in main OSD menu Can be used to activate the Input hot menu when the main OSD menu is not displayed Can be used to Lock/Unlock settings for key activation by holding the key for longer than 30 seconds when the main OSD menu is not displayed
	
Left Key	Move between options and adjust values in the main OSD menu Can be used to activate the Backup Source MSG On/Off hot menu by holding the key for longer than 3 seconds when the main OSD menu is not displayed. Wake up in DPMS mode
	
Right Key	Move between options and adjust values in the main OSD menu. Can be used to activate the Display Mode hot menu by holding the key for longer than 3 seconds when the main OSD menu is not displayed. Wake up in DPMS mode
	
Port A Key	Full Screen Mode Used to change to the backup source input specified in A Used to display the selectable OSD menu for backup source input A by holding the key for more than 5 seconds
	PIP Mode Main input selection function
Port B Key	Full Screen Mode Used to change to the backup source input specified in B Used to display the selectable OSD menu for backup source input B by holding the key for more than 5 seconds
	PIP Mode Sub input selection function

2.2.4 Input/Output Signal Ports

2.2.4.1 PAX-327-C8A-NAS, PAX-327-C8A-NAR

1. DC switch
2. DC power input: +24V power connector
3. DC power output: +5V/+12V power connector
4. RS-232C
5. DVI input/output
6. HDMI input/output
7. DP input
8. USB up

2.2.4.2 PAX-327-C8E-NAS, PAX-327-C8E-NAR (SDI 12G Optional)

1. DC switch
2. DC power input: +24V power connector
3. DC power output: +5V/+12V power connector
4. RS-232C
5. DVI input/output
6. HDMI input/output
7. DP input
8. USB up
9. SDI input/output (optional)

2.2.4.3 PAX-327-C8D-NAS, PAX-327-C8D-NAR (SDI Quad 3G Optional)

1. DC switch
2. DC power input: +24V power connector
3. DC power output: +5V/+12V power connector
4. RS-232C
5. DVI input/output
6. HDMI input/output
7. DP input
8. USB up
9. SDI input/output (optional)

2.2.4.4 PAX-327-C8F-NAS (SDI Quad 3G, Ndcoder Optional)

1. DC switch
2. DC power input: +24V power connector
3. DC power output: +5V/+12V power connector
4. RS-232C
5. DVI input/output
6. HDMI input/output
7. DP input
8. USB up
9. SDI input/output (optional)
10. SFP+ (optional)
11. Status (optional)
12. LAN (optional)

2.2.4.5 PAX-327-C8A-HAS, PAX-327-C8A-HAR, PAX-327-C8A-HAR-AI

1. DC switch
2. DC power input: +24V power connector
3. DC power output: +5V/+12V power connector
4. RS-232C
5. DVI input/output
6. HDMI input/output
7. DP input/output
8. USB up

2.2.4.6 PAX-327-C8E-HAS, PAX-327-C8E-HAR (SDI 12G Optional)

1. DC switch
2. DC power input: +24V power connector
3. DC power output: +5V/+12V power connector
4. RS-232C
5. DVI input/output
6. HDMI input/output
7. DP input/output
8. USB up
9. SDI input/output (optional)

2.2.4.7 PAX-327-C8D-HAS, PAX-327-C8D-HAR (SDI Quad 3G Optional)

1. DC switch
2. DC power input: +24V power connector
3. DC power output: +5V/+12V power connector
4. RS-232C
5. DVI input/output
6. HDMI input/output
7. DP input/output
8. USB up
9. SDI input/output (optional)

2.2.4.8 PAX-327-C8L-HAS, PAX-327-C8L-HAR (SDI Quad 3G, SDI 12G Optional)

1. DC switch
2. DC power input: +24V power connector
3. DC power output: +5V/+12V power connector
4. RS-232C
5. DVI input/output
6. HDMI input/output
7. DP input/output
8. USB up
9. SDI input/output (optional)
10. SDI input/output (optional)

2.2.4.9 PAX-327-C8F-HAS (SDI Quad 3G, Ndcoder Optional)

1. DC switch
2. DC power input: +24V power connector
3. DC power output: +5V/+12V power connector
4. RS-232C
5. DVI input/output
6. HDMI input/output
7. DP input/output
8. USB up
9. SDI input/output (optional)
10. SFP+ (optional)
11. Status (optional)
12. LAN (optional)

2.2.4.10 PAX-327-C8J-HAS (Ndcoder Optional)

1. DC switch
2. DC power input: +24V power connector
3. DC power output: +5V/+12V power connector
4. RS-232C
5. DVI input/output
6. HDMI input/output
7. DP input/output
8. USB up
9. SFP+ (optional)
10. Status (optional)
11. LAN (optional)

2.2.5 Interface on AVAS Back Panel

Name	Function
SFP+	SFP+ port for SFP+ modules only
Status LED	Status LED, will indicate the AVAS operational status: a. Blinking green: operating as normal b. Blinking orange: warning notification c. Continuous orange: FPGA is starting up d. Blinking red: error notification e. Off: no power supplied to the device
1GbE	1GbE (RJ45) has no function yet

Note! *If the power supply cable is unplugged from the power outlet, the AVAS will be disconnected. The power outlet socket should be located near the equipment and be easily accessible.*



2.3 Viewing Angle

The PAX-327-C8 LCD monitor can be vertically adjusted to the preferred viewing angle for maximum comfort.

Caution! *Do not force the monitor past its maximum extension in any direction as this may damage the monitor and the monitor stand.*



2.4 Positioning

Before setting up your workstation, prepare a suitable installation space. This should be a stable, flat, dust-free surface with good surrounding ventilation. Position the monitor screen away from direct sunlight. The glare caused by reflected sunlight may render the screen difficult to read.

Caution!



- When positioning the equipment, ensure that the main ports and sockets are easily accessible.
- Do not place the monitor close to a heat source.
- Do not place the monitor in direct sunlight or near a window. Exposure to direct sunlight or moisture can damage the monitor.


2.5 Connecting AC/DC Power

- Plug the female end of the AC/DC power adapter into the DC power connector.
- Plug the female end of the power cord into the AC power connector on the adapter.
- Plug the male end of the power cord into a power outlet socket.
- The power cable plug may differ according to the standards of each country.


2.6 Connecting Video

- Before connecting the monitor to a PC or other device, ensure all equipment is powered off.
- Connect one end of an HDMI 2.0 cable to the monitor's HDMI port and the other end to the HDMI port of the PC or other device. A DP 1.2 cable can be connected to the monitor's DP port and to the DP port of the PC or other device if equipped with a DP input port.
- Ensure that the cable is securely connected to both the monitor and the PC/ other device. Tighten the connector screws to ensure a secure connection.
- The signal input source can be selected via the OSD menu. Refer to Section 3.1.7 for further details.


PAX-327-C8A-NAS, PAX-327-C8A-NAR




PAX-327-C8E-NAS, PAX-327-C8E-NAR




PAX-327-C8D-NAS, PAX-327-C8D-NAR




PAX-327-C8F-NAS




PAX-327-C8A-HAS, PAX-327-C8A-HAR, PAX-327-C8A-HAR-AI




PAX-327-C8E-HAS, PAX-327-C8E-HAR




PAX-327-C8D-HAS, PAX-327-C8D-HAR





PAX-327-C8L-HAS, PAX-327-C8L-HAR



PAX-327-C8F-HAS



PAX-327-C8J-HAS




Chapter 3

OSD Menu

3.1 OSD Menu Controls


This chapter refers to the On Screen Display (OSD) controls. Using the OSD controls the user can adjust contrast, brightness, display clarity, color temperature and etc. Please read this chapter carefully to get the most out of this monitor.

3.1.1 Display Mode




Menu	Description	Default	Remarks
Display Mode	Used to select one of the various display modes.	1P	

3.1.2 Display Function




Menu	Description	Default	Remarks
Disp Rotate	Used to rotate the screen.	0°	Apply in 1P mode
Disp Latency	Not applicable.		
LR Ratio	Used to adjust the ratio for the left and right displays.		Apply in 2P LR mode
Input Swap	Used to swap input sources for the left and right displays.		
Input Swap	Used to swap input sources for the top and bottom displays.		Apply in 2P TB mode
PIP Position	Used to select one of the various positions on the sub-screen.	Bottom right	Apply in 2P PIP mode
PIP Transparency	Used to adjust the transparency of the sub-screen.	0	
PIP Size	Used to adjust the size of the subscreen.	10	
Input Swap	Used to swap input sources for main and sub-display.		

3.1.3 Select Region




Menu	Description	Default	Remarks
Left Side Right Side Full	Select one region to control: Left side, Right side or Full.	Left side	Apply in 2P LR mode
Top Side Bottom Side Full	Select one region to control: Top Side, Bottom Side or Full.	Top Side	Apply in 2P TB mode
Main Sub Full	Select one region to control: Main, Sub or Full.	Main	Apply in 2P PIP mode
1P 2P 3P 4P Full	Select one region to control: 1P, 2P, 3P, 4P or 4P Full.	1P In	Apply in 4P mode

3.1.4 Picture




Menu	Description	Default	Remarks
Backlight	Used to adjust the backlight luminance.		
Brightness	Used to adjust the screen brightness. Press the < or > button to adjust.	50	Independent of each input
Contrast	Used to adjust the screen contrast. Press the < or > button to adjust.	50	Independent of each input
Sharpness	Used to adjust the screen sharpness. Press the < or > button to adjust.	5	Independent of each input
SDI Format (Option)	Used to select SDI format.	2SI	

3.1.5 Color



Menu	Description	Default	Remarks
Panel Uniformity	Not applicable.		
Gamma	Used to select one of the various gammas.	2.2	Independent of each input
Temperature	Used to select one of the various color temperatures. * Set Gamma to OFF to adjust the red, green, and blue values in the user mode.	Clear	Independent of each input
Color Effect	Used to select one of the various color effects.	Standard	Independent of each input
Profiles	Used to select one of the various profiles. Select the profile and select the Apply icon to load predefined parameters such as Backlight, Input, Brightness, Contrast, Sharpness, Saturation, Hue, Gamma, Temperature, etc. Users can edit each parameter by selecting the Edit icon, and select the Apply icon to load the edited parameter.	OFF	Apply in 1P Mode
Color Format	Used to view the color format of the current input. [RGB, YUV]		
PCM	Not applicable.		
Hue	Used to adjust the screen color hue.	50	Independent of each input
Saturation	Used to adjust the screen color saturation.	50	Independent of each input


3.1.6 Advance



Menu	Description	Default	Remarks
Aspect Ratio	Used to select one of the various screen aspect ratios.	Full	
Over Scan	Used to activate/deactivate the over scan function. [ON,OFF]	OFF	
Switching Time	Used to set the screen switching time.	2	
Over Drive	Not applicable.		
DDCCI	Used to activate/deactivate the DDCCI function. [ON,OFF]	ON	Display Data Channel Command Interface
Ultra Vivid	Used to configure the ultra vivid screen. [OFF, L, M, H]	OFF	
Over Scan Ratio	Used to adjust screen over-scan ratio. * After adjusting the over-scan ratio, set the Over Scan menu from Off to On.	10	
DP Option (Option)	Used to select a DP Version. [D0, D1, D6]	D0	
DP MST (Option)	Used to select a Multi stream transport. [OFF, D0, D1, D6] Clone mode should be set to Off to select DP MST. DP MST is active only if INPUT is selected as DP. Setting INPUT to AUTO SELECT disables DP MST.	OFF	Display Port Multi Stream Transport
DP EDID (Option)	Not applicable		


Clone mode (Option)	Used to activate/deactivate the Clone mode. [OFF, ON] DP MST should be set to OFF to select Clone mode. For Clone mode to operate at 4K_60 Hz, the DP option/D0 on the sub-monitor must also be set to 1.2.	ON	Apply in 1P Mode
Free Sync	Not applicable		
Freeze	Used to activate/deactivate the screen freeze function. [ON,OFF]	OFF	

3.1.7 Input




Menu	Description	Default	Remarks
DP HDMI DVI SDI AVAS Auto Select	Used to select the monitor's display input sources.	Auto	Apply in 1P mode
Left Side	Used to select the input sources for the left and right screens. * See Section 3.1.12 for selectable inputs.		Apply in 2P LR mode
Right Side			
Top Side	Used to select the input sources for the top and bottom screens. * See Section 3.1.12 for selectable inputs.		Apply in 2P TB mode
Bottom Side			
Main	Used to select the input sources for the main and sub screens. *See Section 3.1.12 for selectable inputs.		Apply in 2P PIP mode
Sub			
1P In	Used to select the input sources for 1P, 2P, 3P, and 4P screens. * See Section 3.1.12 for selectable inputs.		Apply in 4P mode
2P In			
3P In			
4P In			

3.1.8 Audio (Option)




Menu	Description	Default	Remarks
Volume	Used to adjust the speaker volume on the monitor.	50	
Mute	Used to activate/deactivate the mute function. [ON, OFF]	OFF	
Stand Alone	Don't Switch Audio Source if Line In is Currently Playing unless a Digital Source is found. (Optional - ready for future use)	OFF	
Audio Source	Used to select the audio source. [Analog, Digital region1]	Digital	

3.1.9 Other



Menu	Description	Default	Remarks
Reset	Used to activate/deactivate the monitor reset function. [YES, NO]	NO	
Menu Time	Used to set the duration of the OSD menu display.	20	
OSD H Position	Used to adjust the horizontal position of the OSD menu.	50	
OSD V Position	Used to adjust the vertical position of the OSD menu.	50	
Language	Used to select OSD language.	English	

3.1.10 Information



Menu	Description	Default	Remarks
Information	Used to access monitor information.		

3.1.11 OSD Input Source

DP HDMI DVI SDI (Optional) AVAS (Optional) Auto Select	Users can select the input source. (Note: The auto select option cannot be selected in multi-display mode.
---	---

3.1.12 PIP Table

Main Source \ Sub Source	DP	HDMI	DVI	SDI (Optional)	AVAS (Optional)
DP	O	O	O	O	O
HDMI	O	O	O	x	x
DVI	O	O	O	O	O
SDI (Optional)	O	x	O	O	x
AVAS (Optional)	O	x	O	x	O

Chapter 4

Safety Information

4.1 Symbols

Symbol	Title
	This symbol is intended to alert the user to the presence of important operating and maintenance (servicing) instructions.
	This symbol is intended to alert the user to the presence of uninsulated “dangerous voltage” of sufficient magnitude to constitute a risk of electric shock.
	This symbol is intended to remind the user to follow the instructions for use, particularly with medical electrical equipment.
	This symbol denotes direct current.
	This symbol denotes the manufacturer.

4.2 FCC Compliance



This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

4.2.1 Class A

This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference. In such cases, users are required to correct the interference at their own expense.

However, any modifications not expressly approved by the manufacturer may void the user's authority to operate the equipment under the FCC Rules.

4.3 CE Marking



This marking is a declaration of conformity in accordance with Article 10(1) of the CE directive.

This product has passed the CE test for environmental specifications. Test conditions for passing include the equipment being operated within an industrial enclosure. To protect the product from damage due to electrostatic discharge (ESD) or electromagnetic interference (EMI) leakage, we strongly recommend the use of CE-compliant industrial enclosure products.

4.4 VCCI (Voluntary Control Council for Interference)



この装置は、情報処理装置等電波障害自主規制協議会（VCCI）の基準に基づくクラス A 情報技術装置です。この装置を家庭環境で使用すると電波妨害を引き起こすことがあります。この場合には使用者が適切な対策を講ずるよう要求されることがあります。

4.5 Disposal of Waste



- Products bearing this mark should not be disposed of with other household waste at the end of their working life. To prevent possible harm to the environment or human health from uncontrolled waste disposal, separate this product from other types of waste and recycle it responsibly to promote the sustainable reuse of material resources.
- Household users should contact either the retailer where they purchased this product, or their local government office, for details of where and how they can take this item for environmentally safe recycling.
- Business users should contact their supplier and check the terms and conditions of the purchase contract. This product should not be mixed with other commercial waste for disposal.

4.6 Safety and Maintenance

Caution! *To ensure optimum performance, follow the safety instructions outlined below when setting up or using the equipment.*



- Do not open the monitor. There are no user-serviceable parts inside. Opening or removing the cover may expose you to dangerous shock hazards or other risks. Refer all servicing to qualified service personnel.
- Do not spill liquid onto the equipment and do not operate the equipment near water.
- Do not insert objects of any kind into the equipment as they may touch dangerous voltage points, which can be harmful or fatal, or may cause electric shock, fire, or equipment failure.
- Do not place heavy objects on the power cord. Damage to the cord may cause shock or fire.
- Do not place the equipment on a sloping or unstable cart, stand, or table as it may fall and cause serious damage.
- To operate the equipment, use an AC power supply cable that complies with international safety standards.
- Do not place any objects on top of the equipment.
- Do not operate the equipment outdoors.
- The fluorescent tube inside the monitor contains mercury. Follow relevant local regulations when disposing of the tube.
- Do not use the equipment in very hot, humid, dusty, or oily environments.
- If the monitor screen is damaged or cracked, handle with care. Do not touch the liquid crystals.
- Ensure adequate ventilation around the equipment to allow heat dissipation. Do not block ventilated openings or place the equipment near a radiator or other heat sources.
- The power cable connector is the primary means of detaching the equipment from the power supply. The monitor should be installed close to a power outlet that is easily accessible.
- Handle with care when transporting. Save the original packaging for moving the equipment between locations.

Caution! *If any of the conditions listed below occur, immediately unplug the monitor from the power outlet and have it serviced by qualified service personnel.*



- The power supply cord or plug is damaged.
- Liquid or objects have penetrated the monitor.
- The monitor has been exposed to rain or water.
- The monitor has been dropped or damaged.
- The monitor does not operate normally according to the operating instructions.

Warning! High Voltage

The AC/DC inverter is equipped with the caution label shown below to indicate the presence of high voltage current. Do not under any circumstances open the monitor to access the inverter. If servicing is required, contact a qualified service technician.



4.6.1 Radio Frequency

- This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to other devices, which can be determined by turning the equipment off and on again, the user is encouraged to try to correct the interference by one or more of the following measures:
 - Reorient or relocate the receiving device.
 - Increase the separation between the equipment.
 - Connect the equipment into an outlet on a circuit different from that to which the other device(s) are connected.
 - Consult the manufacturer or field service technician for assistance.

4.6.2 Keep Away From Windows

- Exposure to rain, moisture, or sunlight can seriously damage the electrical components. To protect the monitor from damage, do not place it near a window.

4.6.3 Storage Temperature

- Do not place the monitor in an environment with a storage temperature of below -20 °C (-4 °F) or above 60 °C (140 °F) as this can cause serious damage.

4.6.4 Image Persistence

- Image persistence is when a residual or ghost image of a previous image remains visible on the screen. Unlike CRT monitors, with LCD monitors image persistence is not permanent, but constant images being displayed for a long period of time should be avoided.
- To alleviate image persistence, turn off the monitor for as long as the previous image was displayed. For example, if an image was on the monitor for one hour and a residual image remains, the monitor should be turned off for one hour to erase the image.

Note! As with all personal display devices, Advantech recommends using a moving screen saver whenever the screen is idle and turning off the monitor when not in use.



4.6.5 Placement and Adjustment of the Monitor

- For optimum performance, a one hour warm-up time is recommended.
- Adjust the monitor height so that the top of the screen is at or slightly below eye level. Your eyes should look slightly downward when viewing the middle of the screen.
- Rest your eyes periodically by focusing on an object farther away and blinking often.
- Position the monitor at a 90° angle to windows and other light sources to minimize glare and reflections. Adjust the tilt of the monitor to avoid reflecting any ceiling lights on the screen.
- Use an anti-glare filter to improve viewing in the presence of reflected light.
- Use a lint-free, non-abrasive cloth for cleaning the surface of the LCD monitor. Do not use any cleaning solution or glass cleaner.
- Adjust the monitor's brightness and contrast setting to enhance readability.
- Avoid displaying fixed patterns on the monitor for long periods of time to avoid image persistence (residual image).

Note! *For optimal viewing, maintaining a 20-inch distance from the screen is recommended.*



Note! *This device generally does not require the user to install software. Any modification of the internal firmware must only be done following specific procedure provided by Advantech.*



Note! *The configuration of the device should mainly be done by the user using the provided hardware keypad.*



Note! *If software tools are provided by Advantech they should be operated using the specific instructions provided. Any software tools provided should only be installed on a computer with appropriate systems to prevent unauthorized access. It is recommended to install anti-virus software and not connect to unsafe external networks.*



Note! *In case of any serious incident that has occurred, please contact the manufacturer and the authorities immediately.*



Appendix A

Cleaning

A.1 Cleaning the LCD Panel

- When the LCD panel becomes dusty or dirty, wipe it gently with a soft cloth.
- Do not rub the LCD panel with cloth of a hard or coarse material.
- Do not apply pressure to the LCD surface.
- Do not use an OA cleaner for cleaning as it may damage or discolor the LCD surface.

A.2 Cleaning the Cabinet

- Ensure the power supply is unplugged and disconnected.
- Gently wipe the cabinet with a soft cloth.
- Dampen a cloth with water and neutral detergent and wipe the cabinet.
- Dry the cabinet using a soft cloth.

Note!



The cabinet surface comprises many plastic materials. Do not clean with benzene, thinner, alkaline detergent, alcohol-based system detergent, glass cleaner, wax, polish cleaner, soap powder, or insecticide. Do not place rubber or vinyl against the cabinet for long periods. These types of fluids and materials can cause the paint to deteriorate, crack, or peel.

Appendix **B**

Troubleshooting

B.1 Non-Responsive Power Button

- Unplug the power cord from the AC outlet to power off and reset the monitor.

B.2 No Picture

- Check that the signal cable is securely connected to the monitor and the display source device.
- Check that the monitor and the display source device are powered on.
- Check that the monitor and display source device are configured to the correct display mode.
- Check that the monitor is compatible with the display card and the settings are all configured correctly.
- Check that the signal cable connector pins are not bent or pushed in.
- Check that the monitor's signal input settings are correctly configured.

B.3 No Video

- If no video is displayed onscreen, turn the monitor power off and on again.
- Ensure that the monitor and display source device are not in power-saving mode.

B.4 Image Persistence

- Image persistence refers to when a residual or "ghost" image of a previous image remains visible on screen. Unlike CRT monitors, with LCD monitors image persistence is not permanent. However, displaying constant images for long periods of time should be avoided.
- To alleviate image persistence, power off the monitor for the same length of time as the image was displayed. For example, if an image was displayed on the monitor for one hour and a residual image remains, the monitor should be powered off for one hour to erase the image.

Note! *As with all personal display devices, Advantech recommends using a moving screen saver whenever the screen is idle and turning off the monitor when not in use.*



B.5 Unstable, Unfocused, or Swimming Image

- Check that the signal cable is securely connected to the monitor and the display source device.
- Check that the monitor is compatible with the display card and the settings are configured correctly.
- If text is presented on screen in a garbled manner, change the display mode and set the refresh rate to 60Hz.

B.6 Incorrect Display Image Size

- Check that the monitor and display source device are configured to the correct display mode (refer to the display card or monitor user manual when changing the graphics mode).



Enabling an Intelligent Planet

www.advantech.com

Please verify specifications before quoting. This guide is intended for reference purposes only.

All product specifications are subject to change without notice.

No part of this publication may be reproduced in any form or by any means, such as electronically, by photocopying, recording, or otherwise, without prior written permission from the publisher.

All brand and product names are trademarks or registered trademarks of their respective companies.

© Advantech Co., Ltd. 2021